

The versatility you've been looking for in arthroscopic positioing.



- Mounts inside or outside the sterile field
- Uses standard rail clamps
- Compatible with all Broadhead Pegboard System pegs and accessories
- Autoclavable and radiolucent polymer provides greater strength and durability
- Complete assembly can be autoclave sterilized
- Made in USA



PATIENT OUTCOMES

Match Grade
Medical
precision surgical equipment

designed by surgeons for surgeons


Soloarc

Field Loc

broadhead

www.matchgrademed.com • (920) 486 -1955

A 3D rendering of a human skeleton, showing the head, neck, torso, and limbs in a surgical position. The skeleton is primarily blue, with the knee and elbow joints highlighted in green, indicating the areas of focus for the surgical equipment.



Most Versatile
Easiest to Clean
Easiest to Use




Uncomplicated
Cleaning

- Entire assembly is completely autoclavable
- Wipeable surface with no hidden crevices that hide organic debris
- Can Mount inside or outside the sterile field

Convenient
Application

- Mounts with standard rail clamps
- Mounts anywhere on the surgical table
- Change position by simply moving a peg
- Adjusts without reaching outside the sterile field



Superior
Versatility

- Kits available for knees or elbows
- Multiple peg lengths and extensions provide an array of positioning options.
- Standard and bariatric accessories available
- Compatible with all Broadhead Pegboard System pegs and accessories