


Stop dealing with pressure ulcers and unnecessary complications.

broadhead



Standard Broadhead System Kit (HP1000)

- Double wide pegs provide 2 more inches of contact area
- Tapered pegboard ends reduce pressure points
- Autoclavable polymer provides greater strength and durability
- Multiple peg lengths provide greater versatility and more secure fit
- Radiolucence polymer allows for comprehensive patient imaging
- Compatible with most standard positioning pegboards
- Made in USA

PATIENT OUTCOMES

Match Grade
precision surgical equipment

designed by surgeons for surgeons

Solorail)

broadhead MINI


field Loc



broadhead



PATIENT OUTCOMES



**Optimal Pressure Relief
Broadest Surface Area
Most Comfortable**

*broad*head



Superior

PRESSURE RELIEF

- Protects skin integrity in lateral decubitus
- Double wide pegs reduce pressure on bony prominences
- Tapered pegboard ends reduce contact with pressure points



Versatile

Utility

- Streamlined assembly as easy as putting a peg in a hole
- Practical alternative for lower extremity arthroscopic thigh bolster and bumps
- Smooth, non-grip surface allows for comprehensive cleaning
- Optional peg shapes and sizes provide support and stability in a broad range of procedures

**Improved
CIRCULATION**

- Broad surface area allows for better blood flow preventing tissue death that leads to pressure ulcers
- Smooth, rounded surfaces can help to avoid shear injury

